Invacare® **Alu Lite**

User guide	EN
Manuel d'utilisation	FR
Manual de Utilizacão	PT
Gebruiksaanwijzing	NL
Manuale d'uso	IT
Manual del usuario	ES
Bedienungsanleitung	DE

*This manual MUST be given to the user of the product.
BEFORE using this product, read this manual and save for future reference*

Yes, you can.

Invacare® **Alu Lite**

User guide	EN	3
Manuel d'utilisation	FR	20
Manual de Utilização	PT	43
Gebruiksaanwijzing	NL	65
Manuale d'uso	IT	87
Manual del usuario	ES	109
Bedienungsanleitung	DE	131

Foreword

The information contained in this manual is subject to change without notice. Some information is submitted under copyright – all rights reserved. Any information in this document cannot be photocopied or duplicated without prior written authorization by Invacare.


As the European and world's leading manufacturer of wheelchairs, Invacare endeavours to supply a wide range of wheelchairs to meet all the needs of the user in everyday life. However, final selection of the wheelchair rests solely with the user and his/her qualified health advisor.

Proper and efficient use of the wheelchair that you have chosen is based upon the medical prescription which was issued for you on the basis of your pathology and the nature of your disability.

Your wheelchair is especially designed to be used indoors, and with certain restrictions outdoor

(Please look at  Warning paragraphs; Part A General Guidelines).

Symbols

Warning!

Indicates a hazardous situation that could result in serious injury or death if it is not avoided.

: Gives useful tips, recommendations and information for efficient, trouble-free use.

 This product complies with Directive 93/42/EEC concerning medical devices.

The launch date of this product is stated in the CE declaration of conformity.

Stamp of the Distributor

Introduction

Dear Customer

Thank you for purchasing an Invacare wheelchair.

This model was designed to provide you with all the benefits and features to meet your needs. Only quality components were selected for your wheelchair based upon rigorous inspections during the entire manufacturing process.

This manual describes the operating limits of your wheelchair, maintenance operations and adjustments that you or your assistant can make.

However, all the repairs (except for inner tubes) as well as some adjustments, require specific technical training and, therefore, must be performed by your distributor.

Intended use

The wheelchair is intended for assisted transport of people who are not able to walk over a long distance.
The wheelchair should only be used on level ground and accessible terrain, as well as indoors.

Warning : Any other or incorrect use of the wheelchair could lead to hazardous situations.

Service life

Even though your Invacare® wheelchair has been designed to provide a long and trouble free life it is inevitable that wear, tear and usage will eventually render the product unusable. Invacare® recommends that the expected service life of this product is five years, providing the product has been correctly maintained according to the manufacturer's recommendations.

TABLE OF CONTENTS

A. GENERAL GUIDELINES

1. Safety and operating limits

- 1.1. Reaching an object from the chair
- 1.2. Transferring to other seats
- 1.3. Tilting
- 1.4. Tilting, kerbs
- 1.5. Stairways

2. Operating instructions

- 2.1. Unfolding and folding the wheelchair
- 2.2. Wheelchair propulsion
- 2.3. Lifting the wheelchair

3. Safety inspection and maintenance

- 3.1. Performance control
- 3.2. General inspection

4. Transportation

5. Warranty

- 5.1. Standard terms and condition
- 5.2. Limitation of liability

6. Summary of operating instructions

B. DESCRIPTION OF YOUR CHAIR

Page		Page
	I. General	
3	1.1. Introduction	13
	1.2. General description	13
	2. Adjustments	
8	2.1. Seat elements	14
	2.1.1 Seat upholstery	14
	2.1.2 Backrest type	14
	2.1.3 Backrest upholstery	15
	2.1.4 Footrest supports	15
9	2.2. Frame	16
	2.2.1 Side frame	16
	2.2.2 Folding system	16
	2.3. Rear wheel	16
	2.3.1 12" wheels	16
	2.3.2 Axles	16
	2.4. Castors	16
	2.5. Brakes	16
	2.5.1 Manual brakes	16
	3 Specification and tool's	
11	3.1. Standard wheelchair specifications	18
	3.2. Tools for adjustments and regular maintenance (Not supplied)	19
	3.3 Dimensions and labelling	19
11		

A. GENERAL GUIDELINES

I. Safety and operating limits

For a safe operation of your wheelchair, the following parameters should be observed :

- Stability and balance


Your wheelchair has been designed to provide the stability you need during normal daily activities.

Any movement in the wheelchair will have an impact on the position of the centre of gravity, which may lead to the wheelchair tipping and a fall.

To improve your safety when you move a lot or you transfer your weight from one place to another, we recommend using seat belts.

- Weight distribution (figure 1)

Many actions cause the user of a wheelchair to reach out, lean over or move about within the wheelchair and outside it. These movements cause a change to normal balance, centre of gravity (G) and weight distribution of the wheelchair.

1

- Weight Limit

The maximum recommended weight of the user is 100 kg. However, the level of activity is an essential factor. For example, an active user who weighs 75 kg may subject the wheelchair to more stress than that of a user who weighs 100 kg. To this purpose, we recommend that you consult your retailer when choosing the model of wheelchair based upon your daily life style.


1.1. Reaching an object from the chair

The limitations on reaching out from a wheelchair indicated in the following diagrams have been calculated based on a representative sample of wheelchair users:

3

- Only the arms should be extended beyond the seat of the wheelchair. (figure 2).

- The body and head should remain within the boundaries of the seat. (figure 3).

2


3


1.1.1 Leaning forward

Do not extend your chest over the armrest (figure 4). In order to reach an object in front of you, you must lean and bend down ; therefore, you must use the castors as a tool (pointing them forward) to maintain stability and balance. An accurate alignment of the wheels is essential for your safety (figure 5).

4

4


5

1.1.2 Leaning backward


Position wheelchair as close as possible to the desired object so that you can simply pick it up by stretching your arm while sitting in the chair in a normal position (figure 6). In any case, do not lean backwards because you may cause the chair to tip (figure 7).

Warning :

The hanging of back packs or similar items onto your chair backposts can affect the rearward stability of your chair, especially when used in combination with recliner backrests. This can cause the chair to tip backwards causing injury. Therefore, Invacare® strongly recommend the use of anti-tippers (available as an option) when using your backposts to carry a back pack or similar items.

6

7

I.2. Transferring to other seats


- If you are more or less able to stand up and if your upper body is sufficiently strong and mobile, you can transfer forward to another seat. Fold the footplate up and push the footrest/legrest to the side, bend your body forward leaning on the two armrests and lift yourself up; then shift your body towards the place where you want to sit while distributing your weight to the arms and hands (figure 8).

As you may find this movement difficult, we recommend that you obtain some assistance from an attendant or carer to ensure the transfer is as smooth and safe as possible.

Lock both brakes and ensure both castor wheels are facing the forward position.

Then, grip firmly the armrest handle, take care of your nails on the skirt-guard edges.

Push on your arms to easily lift yourself up and the attendant can help you to finalise the transfer to the place where you wish to sit.

⚠ Warning :

- Position yourself as close as possible to the place where you wish to sit.
- When transferring, position yourself as far back as possible in the seat to prevent breaking screws, damaging the seat upholstery or causing the wheelchair to tip forward.
- Lock the two brakes.
- Never stand on the footrests when you are getting in or out of the wheelchair (figure 9).

1.3. Tilting (balancing on the rear wheels)

For greater safety, this operation must be performed by an attendant. The attendant should be aware of the required physical effort and use appropriate positioning in order to relieve the strain on his/her back (keep a straight back and bend your knees during this operation).

To tilt the wheelchair, the attendant must firmly grab the handles making sure both are properly fixed. Warn the occupant in the wheelchair before tilting it and remind him / her to lean backwards and make sure that both feet and hands of the user are clear of the wheels.

Place a foot on the footstep tube and apply downward pressure until the chair reaches the equilibrium point. At this stage, the assistant will feel a difference in weight distribution, which usually occurs at approximately 30°. At this point, the wheelchair can get over the obstacle easily.

Finally, the attendant slowly and gradually lowers the front down to the ground, while firmly holding the handles.

Warning :

- Be aware of detachable parts such as armrests or legrests : they must **NEVER** be used as lifting supports as they may be inadvertently released, resulting in possible injury to the user and / or attendant.
- Do not lower the wheelchair suddenly, even if it is several centimetres from the ground, as this may result in injury of the user.

1.4. Tilting, Kerbs



To get on the pavement :

- Method 1 (figure 10)

The attendant positions the wheelchair in front of the pavement facing forward. Attendant tilts the wheelchair backwards until the castors reach the pavement; attendant pushes the wheelchair forward until the rear wheels are against the kerb and again pushes the wheelchair until the rear wheels climb over the kerb.

- Method 2 (figure 11)

In this case, the attendant stays on the pavement and tilts the wheelchair in a backwards position with the rear wheels against the kerb. The attendant tilts the wheelchair backwards until it is balanced and pulls the wheelchair with a steady movement until the rear wheels climb over the kerb ; then he / she lowers the castors, while making sure that the chair is far enough on the pavement so that the castors do not fall into empty space.

To get off the pavement :

The attendant positions the wheelchair facing forward on the pavement and tilts it backwards until it is balanced, then he/she pushes the wheelchair forward until the rear wheels touch the road after getting over the obstacle; then, he / she gradually lowers the castors to the ground.

1.5. Stairways

Because this is a difficult movement, we recommend using two attendants, one in front of the wheelchair and one behind the wheelchair.

To climb stairways (figure 12) :

After tilting the wheelchair to the point of equilibrium, one assistant (at the back) holds the wheelchair up against the first step grasping the handles firmly to lift.

The second assistant, lifts the wheelchair above the stairs, while holding firmly a fixed part of the frame, and holds it while the first assistant takes a step and repeats the operation.

12

The wheelchair must not be lowered until the last step has been passed and the chair is clear of the stairs.

To descend stairways :

Same operation as above, however, reverse steps as above.

Warning :

- Do not attempt to lift the wheelchair by any removable parts (such as armrests, legrests or footrests).
- Avoid using an escalator which may lead to serious injury in the event of a fall.


2. Operating instructions

2.1. Unfolding and Folding the wheelchair

2.1.1 Unfolding the wheelchair (figure 13) :

- With one hand, grab the armrest or the seat support tube on one side of the wheelchair and slightly tilt it towards you (so that the rear wheel and castor lift from the ground) ;
- With the other hand, push on the seat tube until upholstery is fully extended and seat tube sits in receiver.
- Then, engage the two manual brakes, open the footrest/legrest and check the ground clearance (footrest/ground distance - see paragraph B-2.1.5). You can now sit down in the wheelchair.

8

13


2.1.2 Folding the chair (figures 14 & 15) :

- Fold and lock the footrest/legrest toward the front of the wheelchair.
- Swivel the plates into the vertical position. Using both hands, take the centre front and back edges of the seat upholstery and lift it. Or, tilt the wheelchair to one side and close it using the handles on the backrest.

Warning :

- Fold the wheelchair while keeping the seat upholstery upwards to avoid damage by the folding system.
- Caution finger pinch.

14

15

2.2. Wheelchair propulsion


Wheelchair propulsion is provided by attendant only.

2.3. Lifting the wheelchair

First fold the chair (see paragraph 2.1.2), always lift the wheelchair by gripping the frame at the points (A) shown in the figure 16.

⚠ Warning :

Never lift the wheelchair by removable parts (footrests). Ensure the backrest posts are securely locked in place.

16

3. Safety inspection and maintenance

3.1. Performance control

As the user, you will be the first to notice any possible operational defects with your wheelchair. The following table indicates the easiest troubleshooting symptoms to identify and the preliminary inspection that you can perform. In the event that the symptoms persist after adjusting the pressure in the tyres, please consult your retailer. The inner tubes of the wheels are the only components that you can repair yourself (see paragraph B-2.3).

If the screws and nuts are not tight enough or loose, the axle fork angle is not perpendicular to the ground, please consult your retailer to get a competent person to re-set them to the correct position and tighten appropriately. (Recommended torque minimum : M5 = 3-6 Nm; M6 = 7-12 Nm & M8 = 10-20 Nm)

9

The wheelchair swerves to the right	The wheelchair swerves to the left	The wheelchair turns or moves slowly	The castors lift	Creaking and clinking	Play in the wheelchair	Inspections
		●	●	●	●	Make sure that the bolts are tight
●	●	●				Check the fork angle
●	●		●			Make sure that the 2 castors come in contact with the ground at the same time

3.2. General inspection

Your distributor, who has the required technical expertise, is responsible for any wheelchair repairs.

We recommend that you take the wheelchair to your retailer at least once a year for a complete inspection.

Regular maintenance allows the identification and replacement of defective and worn parts, which improves the daily operation of your wheelchair.

Regular inspections to be performed by you or your assistants :

10

a. General

Make sure that the wheelchair folds and unfolds easily.

Make sure that the wheelchair moves in a straight line. (no resistance or deviation)

b. Attendant brakes

Make sure that the attendant brake shoes or mechanism do not touch the moving tyres.

Make sure that the attendant brakes operate easily.

Make sure that the joints are not worn and do not have excessive play.

c. Folding system

Check the folding system for worn or distorted parts.

d. Skirtguard/armrest upholstery

Make sure that all the fittings are properly tightened.

e. Armpad

Make sure that the armpads are in good condition.

f. Seat and backrest upholstery

Make sure that the upholstery is in good condition.

g. Rear wheels

Make sure that the wheel nuts and precision bearings are tight. Make sure that the wheels are parallel to the frame.

h. Castors

Make sure that the axle is tight by turning the wheel the wheel must gradually come to a stop.

i. Fork/steering tube

Make sure that all the fittings are well tightened.

j. Solid tyres

Check the wear of the solid tyre tread.

k. Cleaning and Disinfection

Cleaning: Use only damp clothes and gentle detergent. Do not use abrasive or scouring liquid. Do not use high pressure cleaning devices on ball bearings (front & rear wheels, fork axles)

Disinfection: Spray or wipe disinfection using a tested and recognised product is permitted. A list of the current permitted disinfectants is available from the Robert Koch Institute at <http://www.rki.de>

Make sure you dry the wheelchair if it is wet (e.g. after washing it or going out in the rain).

4. Transportation

The wheelchair is not designed or tested to be used as a seat in a vehicle, it is necessary to use a vehicle seat adapted to your transportation.

When you lift your wheelchair, for example to store it in the trunk of your car, make sure you use a good body posture which doesn't strain your back.

5. Summary of warranty terms

5.1. Standard Invacare terms and conditions

This is to certify that your manual wheelchair is warranted by Invacare for a period of 2 years for the frame, crossbars and all others parts, subject to the following conditions :

- The manufacturer will not accept responsibility for damage caused by misuse or non-observance of the instructions set out in the user manual.
- During the period of warranty, any parts that have become defective due to faulty workmanship or materials, will be renewed or repaired without charge by the Invacare dealer/supplier.
- The warranty will be forfeited should any unauthorised alteration be made to the equipment.
- The Purchaser's statutory rights under the Consumer Protection Act are not affected.

11

5.2. Limitation of liability

This warranty does not extend to the consequential costs resulting from fault clearance, in particular freight and travel costs, loss of earnings, expenses, etc.

Invacare shall not be liable for :

- Natural wear and tear.
- Inappropriate or incorrect use.
- Defective assembly or setting-up by the purchaser or third parties.
- Defective or neglectful treatment
- Use of unsuitable spares.

6. Summary of operating instructions for optimal safety

- Maximum user's recommended weight : 100 kg.
- Do not attempt to reach objects if you have to move forward in the seat.
- Do not attempt to pick up objects from the floor by reaching down between your knees.
- Do not lean over the top of the upholstery back to reach objects located behind you : this may cause you to tip over
- Always engage both attendant brakes simultaneously.
- Do not carry in the stairways or escalator, user seated in the wheelchair whith only one attendant. This may cause serious injury.
- Do not expose the wheelchair to a temperature higher than 40°C.
- To avoid injury, keep your fingers away from mobile parts (folding system, footrests), and maintain good posture before lifting the wheelchair.

⚠ Warning :

Avoid riding on wet areas as well as gravel, grass, etc. (sand and sea water particularly damage ball bearings).

B. DESCRIPTION OF YOUR WHEELCHAIR

I. PRESENTATION

I.1. Introduction

Your wheelchair has been factory set before you purchased it. However, it must be specifically adapted to your needs. The following detailed paragraphs describe the various functions and possible adjustments as well as available options. You can make some adjustments yourself, while others can be made only by your dealer.

Important: based upon the selected model or options, your new Alu Lite wheelchair may be equipped with all of the components or options which are described in the following pages.

Some models or options may be not be available in your country.

 This is a warning symbol, it is essential that the instructions in these paragraphs are followed in order to prevent personal injury to yourself and others around you.

 This is an information symbol about possibility to contact your Dealer for more information.


I.2. General description (see photo)

Your wheelchair is made of various parts and this manual describes only the main parts. We recommend that you become acquainted with the following terms in order to better understand your wheelchair operation :

- **The seat** consists of the seat and backrest upholstery, the backrest and armrests. This unit is designed to provide optimal comfort.
- **The swing-away footrest support** : this is the supporting part between the frame and the footrest which swivels to facilitate transfers and can be removed during transport.
- **The footrest** consists of an adjustable tube and the footplate which supports the foot.
- **The folding frame** consists of side frames and a folding system including the seat rails. These parts constitute the frame, which is the supporting component of the wheelchair and its strength is well tested (checked at 100 kg).
- **The rear wheel** consists of the wheel and axle. The rear wheels ensure the rear stability. They are mounted on fixed wheel support brackets.
- **The castor** consists of the front wheel and the fork. The castors provide front contact with the ground and determine the steering by the direction of the forks.

- **The attendant brake** is a parking brake.

The two attendant brakes are used to secure the wheelchair when stationary.

2. Adjustments

2.1. Seat elements

2.1.1 Seat upholstery

It provides comfortable support to the user.

Standard padded seats are not adjustable; in the event that they become slack, it is recommended that you ask your dealer to replace them.

Keep your fingers away from movable parts to prevent injuries during folding and unfolding !

2.1.2 Backrest types

To prevent falls or possible injury to the user and/or attendant as the backrest canes are used as a support to tilt and/or climb kerb or stairways, make sure that the backrest canes are securely locked in place.

- Folding backrest (photos I and IA)

To save space during transport, operate lever (A) by pulling or pushing it and fold the top of the backrest.

To return to the initial position, bring the top in the vertical position; it locks automatically.

⚠ Always make sure that the backrest is properly locked in place before the user settles down in the wheelchair to prevent any injuries !

2.1.3 Backrest upholstery


- The padded backrest provides excellent daily comfort for the user who does not need specific support for the upper body.

⚠ In the event that the upholstered backrest slackens, ask for a replacement from your Dealer.

2.1.4 Footrest supports

- Standard footrest supports (photo 2) swing away during transfers and can be removed during transport.

Operate lever (A) by pushing sideways and swivelling towards outside. To return to the initial position, align the footrest support it locks automatically.

To remove the footrest support, simply pull up after unlocking the assembly.

To refit it, align the holes on the side mounting with the support and press on, while still in the unlocked position



- Footrests (photo 3) : the footplate can be lifted during transfers footrests are height adjustable. Loosen the screw (B) to adjust to the desired height, firmly tighten the screw after adjustment.

Straps for standard legrest : to ensure a good position of the feet.

Note : the standard footrest supports and legrests are mounted in pairs on the wheelchair; whenever you remove them, remember that you have a right side and a left side !

⚠ Never lift the wheelchair by the footrest supports or legrests !

Keep your fingers away from movable parts during folding, disassembling or adjustment to prevent injuries !

2.2. Frame

2.2.1 Side frame

The side frames are designed to accommodate fixation of the front and rear wheels. No adjustment is required on side frames.

2.2.2 Folding system

It consists of two cross-bars which integrate the seat rails.
To fold and unfold your wheelchair, see chapter A “General paragraph 2.1. No adjustment is required on folding system.

2.3. Rear wheels

2.3.1 12" Wheels

The 12" (305 mm) rear wheels are composite wheels. They are delivered with solid tyre.

2.3.2 Axles

The wheelchair is equipped with fixed axles on rear wheels.

2.4. Castors

The solid front wheels are available in 8" x 1 1/4" (200mm x 28 mm) diameter.

2.5. Brakes

2.5.1 Manual brakes

- Locking Hand Brake (photo 4)

1. Push down on the bottom portion of the brake handle (A) until an audible click is heard.

2. Refer to **Unlocking Hand Brake**, see following paragraph.

- Unlocking Hand Brake (photo 4)

1. Pull up on the top portion of the brake handle (A) to release.


2. Release the brake handle.

- Using Hand Brake (photo 4)

1. Pull up on the top portion of the brake handles (A) toward the push handles.

2. Do one of the following:

- Remain stationary - Hold the brake handle up.
- Continue mobility - Release the brake handle.

Adjusting brake handle (photo 5)

Test the brake handle. Observe how the brake lever engages the wheel. If the tension of the brake is too loose or too tight, adjust using the following steps:

1. Loosen the brake handle adjustment nut (A) by turning clockwise.
2. Perform one of the following:
 - Loosen brake handle tension - Turn the brake handle nut (B) clockwise.
 - Tighten brake handle tension - Turn the brake handle nut (B) counterclockwise.

NOTE: If the brake lever is scraping the rear wheel, the brake handle must be loosened.


- Tighten brake handle tension - Turn the brake handle nut (B) counterclockwise.

NOTE: If the brake lever does not respond quickly, tighten brakes.

3. Turn the brake handle adjustment nut (A) counterclockwise to secure in place.

4. Do one of the following:

- Acceptable tension - Repeat STEPS 1 - 3 for the opposite side, if necessary.
- Unacceptable tension - Refer to Adjusting Cable Adjuster see following paragraph.

5

5A

Adjusting Cable adjuster unit (photo 6)

Test the brake with the brake handle. Observe how the brake lever engages the wheel. If the tension of the brake handle is still too loose or too tight, adjust using the following steps:

1. Ensure brake handle is not locked.
2. Do one of the following:
 - To Loosen: Brace the jam nut (C) while turning the adjustment screw (D) clockwise.
 - To Tighten: Brace the jam nut (C) while turning the adjustment screw (D) counterclockwise.
3. Repeat with other brake.


Adjusting Wheel brake (photo 7)

1. Loosen the bolt and locknut that secure the wheel brake assembly to the transport chair frame.


2. Adjust the position of wheel brake until the measurement X between the rear wheel and the wheel brake shoe is between 3 and 4 mm.

3. Securely tighten the bolt and locknut.

4. Engage the wheel brake.

6

5. Push against the transport chair and determine if the wheel brake engages the wheel brake shoe enough to hold the transport chair.
6. Repeat the above procedures until the wheel brake holds the transport chair.
7. Repeat STEPS 1-6 for the opposite wheel brake.

3. Specifications and tool's

3.1. Standard wheelchair specifications

Maximum user weight :	100 kg
Seat width :	40,5 / 45,5 cm
Seat depth :	43 cm
Floor/seat height :	49 cm
Rear wheel :	12" (305 mm) solid tyre
Castors :	8" (200 mm) solid tyre
Parking brake :	Attendant brake
Backrest :	Folding
Armrests :	Fixed
Footrest supports :	Removable and swing-away
Seat upholstery :	Black nylon on reinforced upholstery
Frame :	Aluminium, epoxy coated
Wheelchair average weight :	12 kg

3.2. Tools for adjustments and regular maintenance (not supplied)

Function	Tool
Brake	8 mm open-end Wrench
Footrest tube	13 mm open-end Wrench
Armpad	Screwdriver Pozidriv n°2
Castor	5 mm Allen key 13 mm open-end Wrench

After sale and disposal recommendations

- It is compulsory to use original Invacare® spare parts which you can buy through any Invacare® dealer.
- For repair, please contact your local Invacare® dealer.
- Disposal : the metal parts can be disposed of for scrap metal through recycling. Plastic parts are disposed of as plastic scrap. Disposal must be carried out in accordance with the respective national regulation. Please apply to your municipal authorities/local government for details about local disposal companies.

3.3. Dimensions and labelling

Picture	Description	Min/Max value
	Seat effective width (mm)	405/455
	Overall width (mm)	520/570
	Width of folded wheelchair (mm)	250
	Total height (mm)	945
	Height from ground to back seat (mm)	440
	Height from ground to front seat (mm)	490
	Backrest height (mm)	450
	Wheelchair height when backrest is folded (mm)	750
	Backrest height including headrest (mm)	N/A
	Overall length (mm)	1030
	Length without footrest (mm)	750
	Distance between front wheel and rear wheel (mm)	375
	Backrest angle (0°)	0
	Bracket angle (0°)	60
	Distance between footrest and seat (mm)	340/460
	Distance between armrest and backrest (mm)	355


Picture	Description	Min/Max value
	Fixed armrest height (mm)	230
	Adjustable armrest height (mm)	N/A
	Notched armrest height (mm)	N/A
	Rear wheel diameter (mm)	310 mm
	Handrim diameter (mm)	N/A
	Seat plane angle (0°)	3
	Overall weight when empty (kg)	12
	User's maximum weight (kg)	100
	Turning circle (mm)	890
	Cushion depth (mm)	N/A
	Weight of the heaviest removed part (kg)	10,6
	Seat material Fire class	M4 EN 1021 1-2

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A symbol consisting of an exclamation mark inside a circle.

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The CE mark, which consists of the letters 'CE' enclosed in a circle.

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Invacare® Poirier SAS
Invacare® Ireland Ltd
Invacare® AB
Invacare® B.V.
Invacare® Ltd
Invacare ® France Operations SAS